

martin



Electrosurgery


KLS Martin ME 411

Utmost Safety and Reliability
Thanks to an Advanced Safety Concept

www.klsmartin.com

Progressively increasing/ decreasing power output

In the lower output range, the user can select the desired HF output power value in 1-watt steps, which allows very precise dosing. As the power adjuster is rotated clockwise, output power increases on a non-linear, progressive basis. This means that high power reserves are available in the upper power range. Thanks to these unique power-range



Integrated dynamic, self-adaptive HF power output control


With the help of the dynamic, self-adjusting power output control that has been integrated into the system, the KLS Martin ME 411 automatically selects the most suitable operating point. This guarantees optimum cutting results irrespective of tissue impedance or cutting speeds.



Highest patient safety thanks to Patient Control System (PCS)


Using a dual-pad neutral electrode, the integrated KLS Martin Patient Control System (PCS) continually monitors the neutral electrode for proper application. So if the electrode is only insufficiently in contact with the patient's skin, PCS inhibits HF-power and gives alarm, before critical situations arise.

Moreover, the system also monitors single-pad neutral electrodes in acc. with IEC 601-2-2.













Better protection against mal- operation by PIN Error Control

The KLS Martin ME 411 performs a power-on-self-test routine whenever the system has been switched on. If an error or malfunction is found, HF power output is blocked and the respective error code is displayed.




10 different types of current for a broad range of indications

- 
Monopolar cutting 1
 Pure cutting current allowing smooth cuts, with optional integration of a starting cut pulse.
- 
Cutting URO 2
 Cutting current blended with coagulation current, for use in liquid media (e.g. papillectomy).
- 
Spray coagulation
 Non-contact surface coagulation (fulguration). In tissue contact under liquid media (e.g. TUR) useful for staunching micro-bleedings.
- 
Monopolar cutting 2
 Cutting current blended with coagulation current.
- 
Contact coagulation 1
 Coagulation current with very deep-reaching effects; electrode in direct contact with the tissue.
- 
Bipolar cutting 1
 Pure bipolar cutting current, to be used with sling-type electrodes, for example.
- 
Cutting URO 1
 Cutting current allowing smooth cuts in liquid media (e.g. TUR), with optional integration of a starting cut pulse.
- 
Contact coagulation 2
 Coagulation current with deep-reaching effects; electrode in direct contact with the tissue. Particularly suitable in TUR.
- 
Bipolar cutting 2
 Bipolar cutting current blended with coagulation current to be used f. e. with electro-surgical forceps, for example.
- 
Bipolar coagulation
 Automatic bipolar coagulation to be used with bipolar tweezers or bipolar instruments.




Highest safety for the patient by patented active leakage current compensation

The KLS Martin ME 411 monitors and actively controls the HF leakage currents and compensates them. This helps to avoid the risks associated with so-called „accidental HF burns“, thus providing additional safety to the patient.



Advanced safety concept based on two microprocessors working in parallel

The tasks of monitoring and controlling the various system components are shared by two microprocessors working in tandem. As these two modules also monitor each other, they provide the basis for KLS Martin's highly advanced safety concept.



Bipolar Coagulation with AUTO-START and AUTO-STOP Function

The bipolar coagulation current is activated automatically by tissue contact and will automatically cut off the bipolar HF current as


soon as the coagulation process has been completed, prior to unwanted side-effects like tissue glueing or pollution.



Multi-function connector for connecting intelligent accessories

The bipolar output port of the KLS Martin ME 411 has been designed as a multi-functional interface for connecting intelligent accessory equipment. Thus, this interface provides


for proper information transfer between the HF generator and the connected accessory part. Being upward compatible, this interface also takes account of future developments.



Electrosurgical Unit ME 411

Utmost Safety and Reliability

Thanks to an Advanced Safety Concept



Advantages of the KLS Martin ME 411

- Utmost patient safety by constant monitoring and active compensation of HF-leakage currents
- Constant monitoring of neutral electrode for proper application via PCS
- Safety of operation due to power-on-self-test, including detailed error messages
- Digital setpoint value display in accordance with the selected output power level
- Integrated dynamic, self-adaptive HF output control for reproducible cuts
- Progressive power output adjustment
- 10 different types of current available for monopolar and bipolar application
- Two independent power generators
- Two parallel μ -processors
- Two URO-CUT Functions particularly for underwater-cutting
- Serial interface for MABS
- Comfortable bipolar Cut or Coagulation
- Bipolar Autostart/Autostop function
- Symmetrical Bipolar Output, significantly reducing tissue glueing or pollution of electrodes
- Bipolar instrument recognition for use of intelligent accessories
- Broad range of high-quality, practice-oriented accessories – see Martin Accessories catalog



KLS Martin Argon-Beamer-System (MABS)

In combination with the KLS Martin Argon Beamer MB 181 the KLS Martin HF-unit ME 411 with its intelligent interactive interface forms the KLS Martin Argon-Beamer system. The MABS inaugurates new applications for superficial soft coagulation of tissue in open or endoscopic surgery.



KLS Martin ME 411 Specifications

Supply voltage	100 V / 115 V / 127 V / 230 V / 240 V, 50 Hz - 60 Hz Adjustable by the Technical Service		
Power input	in standby mode: 1 VA with no HF power output: 42 VA at max. power output: 880 VA		
Class of protection	I		
Classified acc. to MDD	IIb		
Leakage currents	in acc. with IEC 601, Part 2-2		
Type of equipment	CF; defibrillator-proof		
Nominal frequency	400 kHz		
Pulse frequencies	30 kHz at cut 2, URO cut 2 and contact coagulation 1 65 kHz at contact coagulation 2 50 kHz at spray coagulation 30 kHz at bipolar cutting 2		
HF Power		CF	
Cutting 1	max. 320 W at 350 ohms	1.8	max. 2700 Vss
Cutting 2a	max. 320 W at 350 ohms	2.3	max. 3500 Vss
Cutting 2b	max. 320 W at 800 ohms	2.6	max. 4300 Vss
Urolog. cutting 1	max. 320 W at 350 ohms	1.8	max. 2700 Vss
Urolog. cutting 2	max. 320 W at 800 ohms	2.6	max. 4800 Vss
Contact coagulation 1a	max. 250 W at 200 ohms	1.8	max. 1300 Vss
Contact coagulation 1b	max. 250 W at 200 ohms	3.0	max. 2600 Vss
Contact coagulation 2	max. 150 W at 500 ohms	5.4	max. 4300 Vss
Spray coagulation	max. 120 W at 1000 ohms	5.5	max. 6000 Vss
Bipol. cutting 1	max. 80 W at 500 ohms	1.8	max. 900 Vss
Bipol. cutting 2	max. 80 W at 500 ohms	2.1	max. 1200 Vss
Bipol. coagulation	max. 80 W at 100 ohms	1.8	max. 400 Vss
Duty type	intermittent (10s/30s = duty factor of 25%)		
Sound level	HF activation 55 dB(A), (adjustable 50 dB - 60 dB by Technical Service) Alarm: 65 dB(A)		
Weight	14.2 kg		
Interference suppression	limits in conformity with EN 55011 interference immunity in conformity with IEC 801 CE marking conform with 93/42/EEC		
Dimensions	405 mm x 135 mm x 380 mm (WxHxD)		

Ordering Data

80-041-01	KLS Martin electrosurgical unit ME 411, with main cable, no accessories included
80-140-00	Manual-switch accessories set, large
80-140-01	Foot-switch accessories set, large
80-140-02	Manual-switch accessories set, small
80-150-00	Accessory Set bipolar

martin

Gebrüder Martin GmbH & Co. KG

A company of the KLS Martin Group
Ludwigstaler Str. 132 · D-78532 Tuttlingen
Postfach 60 · D-78501 Tuttlingen
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com

04.07 · 90-657-02 · Printed in Germany
Copyright by Gebrüder Martin GmbH & Co. KG
Alle Rechte vorbehalten.
Technische Änderungen vorbehalten.
We reserve the right to make alterations.
Cambios técnicos reservados.
Sous réserve de modifications techniques.
Ci riserviamo il diritto di modifiche tecniche.