

BEACH CHAIR



Scandia Back plate module for shoulder surgery, electrically adjusted with removable side plates. 3D adjustable head support with comfort cushions. Code 60280




Side support for Shoulder surgery device, code: 602801



Straps for head support (12pcs package). Code 602804

OPERATING TABLE SECTIONS



# SCANDIA<sup>TM</sup>

## OPERATING TABLE

### ACCESSORIES

# SCANDIA<sup>TM</sup>

## OPERATING TABLE

### ACCESSORIES

This competitive selection of accessories is designed to complement operating tables in surgical procedures. These high-quality accessories ensure the tables to be suitable for many specialized procedures. All materials and components used have been carefully selected and manufactured. The safe and easy fixing and ergonomic

design provide efficient and time saving surgical process. The well formed surfaces provide comfort and safety for the patient. Most of Lojer Accessories are made for standard EU side rails (10 x 25 mm), so they are suitable for any Operating Table Brand.

#### TABLE MATTRESSES



Scandia Comfort mattresses have a unique double layer structure, in which the upper visco-elastic layer forms according to patient's forms and body temperature. Velcro – attachment.

Mattress, head section 75 mm	Code 14843
Mattress, head section 120 mm	Code 14843120
Mattress, 25-section, 75 mm	Code 14461
Mattress, 25-section, 120 mm	Code 14461120
Mattress, 25-section with gynaecologic cutout, 75 mm	Code 14463
Mattress, 25-section with gynaecologic cutout, 120 mm	Code 14463120
Mattress, 40-section, 75 mm	Code 14460
Mattress, 40-section, 120 mm	Code 14460120
Mattress, 40-section with gynaecologic cutout, 75 mm	Code 60155
Mattress, 40-section with gynaecologic cutout, 120 mm	Code 60155120
Mattress set with 3 sections, 75 mm	Code 60180-2
Mattress set with 3 sections, 120 mm	Code 60182-2
Mattress set with 6 sections, 75 mm	Code 60180
Mattress set with 6 sections, 120 mm	Code 60182
Long mattress with separate head section 75 mm	Code 60185
Long mattress with separate head section 120 mm	Code 60187

#### ANESTHESIA SCREENS AND INFUSION STANDS



Anaesthesia frame, clamp included, code 10310. Optional extension arm available, code 60328



Anaesthesia frame, flexible, clamps included. Code 10311



Infusion stand Ø 18 mm, stainless steel. One hand height adjustment. Code: 60120. Mobile Infusion stand, stainless steel, castors 5 x 80 mm. One hand height adjustment. Code: 60130



#### CLAMPS



Rotary rod clamp, stainless steel. Code: 61283



Rapid rod clamp, aluminium. Code: 61264



Easy lock blade clamp for accessories with flat mounting blades. Code: 60200

#### ARM RESTS AND TABLES



Armboard with one hand balljoint release, clamp included, without cushion. Code 10380. Model 6010387 includes also a rotatable clamp function.



X-ray transparent Armboard, clamp included. Code: 10385



Comfort cushion (concave) for armboards, code 10382H



Concave cushion for armboards, injection molded. Code 10382



Arm and hand surgery table with detachable supporting leg, clamps included. Code 10390



Multi-Task Armboard, lockable ball socket combined with 15 cm of horizontal travel. Code 60322.

#### SUPPORTS AND STRAPS



Table width extender (pair, includes clamps and cushions), length 500mm. Code 10423



Lateral support, wide. Code 51461416



Lateral support, narrow. Code 61406



Lateral support, circular swivelling soft pad. Code 603211. Requires clamp 60200.



Lateral support, rectangular swivelling soft pad. Code 603212. Requires clamp 60200.



Supporting Roll, vertical. Code: 61421



Supporting Roll, horizontal. Code 10369

SUPPORTS AND STRAPS



Shoulder Support, pair. Code 61253



Foot supports, pair. Code 7510175



Legstrap with clamps, pair, clamps included. Code 10372



Padded Leg and Body strap, wide. Code 10464



Patient restraint strap with Velcro. Code R7842



Wristclamp with strap, clamp included. Code 10370



HEAD RESTS



Flexible head support for ENT-, plastic and eye surgery. Code 60327. Requires adapter 60110



Adapter for special head rests, Mayfield® compatible. Code 60110



Circular head cushion. Code 10333



Head rest for neurosurgery. Code 6010335. Requires adapter 60110

OTHER ACCESSORIES



Laminectomy frame, pair of cushions with abdominal cut out, manually operated, code: 6010485



Uro Catcher System (code 60324) provides excellent fluid containment. Flexible drape support frame provides unobstructed surgeon access during manipulation and springs back to its original state. Compatible with all suction systems for accurate fluid measurement and metering. Sterile bags (10 pcs), code 60235.



Perineal post with cushion. Code 60115. Requires clamps 61264.



ORTHOPAEDIC ATTACHMENT



Orthopaedic extension device for orthopaedic operations and traumatological treatment of fractures.  
Code 60290



Leg board with cushion. Code 61120110



Clamp for orthopaedic attachment.  
Code ES106



Tilting tractor for orthopaedic attachment.  
Code 60291



Extension bar accessory rail. Code ES15334



Scandia Foot control unit, three programmable functions. Code 60270

LEG RESTS AND SECTIONS



Gas spring assisted 3D soft pad stirrups, pair. stirrups allow for easy adjustment of abduction and lithotomy while maintaining the sterile field. They enable safe and easy positioning while providing enhanced surgical site access. Code 60320. Requires clamps 60200



Leg holders with ball joint, pair. Code 10450



Padded strap for leg holders. Code 10452



Arthroscopic leg holder. Code 6128920



Split Leg Section, pair. Code 60160



Scandia Gynecologic/urologic section 400 mm (female/female). Code 60150



Divided leg section, dual joint abduction. Code 60161