




STARKSTROM

patient equipment transfer trolley


Quality Medical Equipment and Service



230VAC BUT-GOING FEED TO FIN... SOCKETS
2 x 4mm... MINIMUM
230VAC INC... LY
2 x 10mm... RUM
COMMON ALARM TO ALARM UN... (s)
3 x 0.75...
COMMON ALARM (vf) FROM...
2 x 0.75mm... SCREENED



Excellence & Innovation in Medical Engineering



transfer trolley






The Starkstrom Patient Equipment Transfer Trolley allows for intra-hospital patient transfers without disruption in care and improves patient safety while reducing valuable transfer setup time. All patient life support equipment can be securely attached to the trolley thus avoiding the expensive equipment damage that commonly occurs when equipment is transported on the patient bed.



Unique universal
bed docking
mechanism



dimensions



key benefits

- Supports all essential equipment necessary for an inter-departmental transfer
- A single "plug-in" point allows permanently mounted equipment to be re-charged
- Quickly and easily docks with all beds to provide a single manageable unit
- Optimum manoeuvrability for patient bed and docked trolley in the tightest corridors
- Eye level positioning of equipment increases ease of patient monitoring whilst in transit
- Rapidly reduces time required for any inter-departmental transfer
- Anti-static wheels and RCD prevents electrical shocks

technical specifications

- Medical Gas Cylinder saddles suitable for sizes A-F
- Electrical outlet strip for plugging in life support equipment
- 5 Metre automatic reel fitted to an RCD protected 13A fused plug
- 4 Anti-static swivel castors with brake lever fitted to rear wheels
- 2 Consumable drawer unit
- 3 Integral utility posts
- 1 Mid level shelf with securing straps
- 2 Bed docking attachments
- 2 Manoeuvring handles
- 1 Vertical utility poles with drip attachment to one
- 1 Central utility pole (vertical)
- 1 Height adjustable infusion pole with a 4-hook drip stand
- Bracket to hold Oxylog Ventilator ±
- Oxylog Ventilator Mounting Plate ±
- GCX Monitor Arm ±
- BREAS LVT 1000 Respirator Holster ±
- Neurological MPM -1 Monitor Frame ±
- TRAM-RAC Module Cage Mounts ±
- Maquet SERVO i Ventilator Bracket ±




± **Optional Extras. Manufactured to individual specification to suit existing hospital equipment.**



British Company

- Clinical Lighting / Cameras and Battery Backups
- Media Solutions / Audio Visual Systems / Cameras and Monitors
- Clinical Medical Gas Pendants / Beams
- Operating Tables
- Interdepartmental Patient Transfer Units
- Surgeons Theatre Control Panels (Touchscreen / Membrane / Stainless Steel)
- PACS Viewing Stations
- Isolated Power Supplies (IPS)
- Uninterruptible Power Supplies (UPS)
- Ultra Clean Ventilation Canopies / Laminar Flow
- Equipotential MEIGaN Compliant Earth Reference Bars / Equipotential Earth Sockets / Clean Earth Blue Sockets
- Illuminated Warning Signs
- Medical Consumables
- Time of Day / Elapsed Time Clock Systems
- Maintenance and Support



STARKSTROM

256 Field End Road Eastcote
Ruislip Middlesex HA4 9UW

T: +44 (0) 208 868 3732

F: +44 (0) 208 868 3736

E: info@starkstrom.com

www.starkstrom.com