



*mar*Test

*... makes sure that only properly functioning
HF accessories find their way into the operating theater.*

- Safe
- Versatile
- Easy-to-use



- Safe
- Versatile
- Easy-to-use

The test unit "marTest" was developed for testing HF-surgical accessories (i.e. monopolar/bipolar cables, electrode handles, neutral electrode cables, KLS Martin IQ instruments).

The unit allows for extremely easy functional testing. The unit is equipped with connections for electro-surgical equipment made by various manufacturers. This allows testing HF-surgical accessories from various manufacturers. Faulty HF accessories can therefore be eliminated in advance.

This ensures that faulty accessories do not find their way into the operating theater in the first place.

Ordering information

80-105-00-04
Test unit marTest

Consisting of:

- Test unit marTest
- Power unit with European primary connector
- Set of cables for insulation testing for instruments with 4-mm socket
- Storage case

80-105-01-04
Cable for insulation testing for bipolar instruments

Gebrüder Martin GmbH & Co. KG
A company of the KLS Martin Group
KLS Martin Platz 1 · 78532 Tuttlingen · Germany
Postfach 60 · 78501 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com