




marLED[®] E1

Examination light



Reliable and safe –
designed to simplify the day-to-day
work of specialists

Reliable surgical lighting is indispensable in clinics and doctors' surgeries. The examination lights from KLS Martin make day-to-day work easier both for examinations and surgery preparations.



Efficient LED technology and flexibility
for optimum support
in day-to-day work

These examination lights bring state-of-the-art LED technology to the field of medical examinations and stand out with high flexibility and good positioning of the light head in day-to-day work environments.

Variable mounting options –
adaptable for any room

Ceiling mounting, wall mounting, rail clamp or mobile tripod: marLED® E1 examination lights can be used flexibly and will cover all of your room lighting requirements.





Many advantages
and benefits for you

- State-of-the art LED technology
- Exceptional lighting power and efficiency
- Long LED service life
- Ease of handling
- Functional and ergonomic design
- Comfortable positioning
- Ergonomic handle
- Low weight
- Completely enclosed light system
- Easy to clean
- High standards of hygiene

Advantages for your
ecological balance

For everyone who is already thinking about tomorrow today, marLED® E1 is the lighting solution of choice. Not only because it ensures optimized lighting of the wound area, but because it also delivers impressive energy efficiency. This certainly adds up. For the clinic, for patients – and for the environment.



marLED® E1 models



Mobile light



Wall light

Ceiling light

Individually tailored to your needs

Whether as a ceiling light, wall light or mobile light, or with an adapter for standard rails, marLED® E1 examination lights help to make work processes as smooth as possible for doctors and care staff, and they can be adapted extremely well to any conditions in the room.

Ordering data

Item no.	marLED® E1 examination light
89-701-00-04	marLED® E1 mobile
89-701-01-04	marLED® E1 wall
89-701-02-04	marLED® E1 ceiling
89-701-03-04	marLED® E1 rail attachment

Technical data

	marLED® E1
Central illuminance (0.5 m distance)	50,000 Lux
Color temperature	4,400 K
Color rendering index [Ra (1-8)]	> 93
Illumination field diameter d10 (0.5 m distance)	17 cm
Supply voltage	100-240 V; 50/60 Hz
Power rating	20 W

For production reasons and metrological reasons, all of the data provided relating to the light output has a tolerance of +/- 5%.

KLS Martin Group

KLS Martin France SARL
68200 Mulhouse · France
Tel. +33 3 89 51 3150
france@klsmartin.com

KLS Martin UK Ltd.
Reading RG1 3EU · United Kingdom
Tel. +44 1189 000 570
uk@klsmartin.com

KLS Martin do Brasil Ltda.
CEP 04.531-011 São Paulo · Brazil
Tel.: +55 11 3554 2299
brazil@klsmartin.com

Gebrüder Martin GmbH & Co. KG
Representative Office
121471 Moscow · Russia
Tel. +7 499 792-76-19
russia@klsmartin.com

Martin Italia S.r.l.
20864 Agrate Brianza (MB) · Italy
Tel. +39 039 605 6731
italia@klsmartin.com

Nippon Martin K.K.
Tokyo 113-0024 · Japan
Tel. +81 3 3814 1431
nippon@klsmartin.com

KLS Martin Australia Pty Limited
Artarmon NSW 2064 · Australia
Tel.: +61 2 9439 5316
australia@klsmartin.com

Gebrüder Martin GmbH & Co. KG
Representative Office
201203 Shanghai · China
Tel. +86 21 5820 6251
china@klsmartin.com

Martin Nederland/Marned B.V.
1271 AG Huizen · The Netherlands
Tel. +31 35 523 4538
nederland@klsmartin.com

KLS Martin L.P.
Jacksonville, FL 32246 · USA
Tel. +1 904 641 7746
usa@klsmartin.com

KLS Martin SE Asia Sdn. Bhd.
14100 Simpang Ampat, Penang · Malaysia
Tel.: +604 505 7838
malaysia@klsmartin.com

Gebrüder Martin GmbH & Co. KG
Representative Office
Dubai · United Arab Emirates
Tel. +971 4 454 16 55
middleeast@klsmartin.com

Gebrüder Martin GmbH & Co. KG
A company of the KLS Martin Group
KLS Martin Platz 1 · 78532 Tuttlingen · Germany
Postfach 60 · 78501 Tuttlingen · Germany
Tel. +49 7461 706-0 · Fax +49 7461 706-193
info@klsmartin.com · www.klsmartin.com

